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[English tutorial](#)

[Instructions en français](#)

Remarque importante:

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The screenshot shows the SpecQuest web interface. Key elements include:
- User information: Group (OC1), Name (Person, Muster), Phone (1234), e-Mail (muster_person@uni-bielefeld.de), Expiration (31.08.2021).
- Submission Form: Mass Spec Submission Form (selected), NMR Submission Form (unselected).
- Sample Data: Sample Name (MP-01), Purification (- none - / mixture / raw), melts/boils (don't know).
- Empirical Formula: C8H10N4O2.
- Toxicity (4): harmful (checked), carcinogenic, explosive, irritant, toxic, untested.
- Properties (4): light sensitive, moisture sensitive, oxidation sensitive, small amount (< 100µg), stable (checked), strong odor.
- Sample is soluble in: (4): Acetone (Acetone), Acetonitrile (ACN), Chloroform (CHCl3), Cyclohexane (CyH), Dichloromethane (Cl2), Ethanol (EtOH).
- Analytical Requirements: Request for: (2): Absence of Compound, Accurate Mass (+5 ppm), Fragmentation Pattern, MW (+/- 0.2 Da) (checked).
- Analogous samples: measured with: (2): CI (Chemical Ionization), EI (Electron Ionization), ESI (Electrospray Ionization) (checked), MALDI (Matr. Ass. Laser Des. Ion.).
- Chromatography necessary: (1): GC-MS, HPLC-MS, TLC Extraction.
- Comment/Info for the Operator: pharmakologisch aktiv!
- Small Reaction Scheme Preview: A chemical structure of a nucleoside.
- Create MS Sample Submission Form button.
- Do not resize in drawing software before insert! (check quality with preview!)
- Insert Reaction Scheme from Clipboard button.
- Clear button.

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Sample dissolved
 Empirical Formula is unknown! a)
Approx. Molecular Weight
300-350 Da b)

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carcinogenic
 explosive
 harmful
 irritant
 toxic
 untested

oxidation sensitive
 small amount (<100r
 stable
 strong odor
 temperature sensitiv

Dichloromethane (Cl
 Ethanol (EtOH)
 Ethylacetate (EtOAc)
 Isopropanol (IsoProp
 Methanol (MeOH)
 Tetrahydrofuran (TH

Analytical Requirements

Request for: (2) Analogous samples have successfully been measured with: (2) Chromatography necessary (1)

Absence of Compound
 Accurate Mass (+5 ppm)
 Fragmentation Pattern
 MW (+0.2 Da)

CI (Chemical Ionization)
 EI (Electron Ionization)
 ESI (Electrospray Ionization)
 MALDI (Matr. Ass. Laser Des. Ion.)

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Routine ESI: suited for polar, high boiling compounds, even salts

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The screenshot shows the 'SpecQuest' software interface. At the top, there are fields for 'Group' (OC1), 'Name' (Person, Muster), 'Phone' (1234), 'e-Mail' (muster.person@uni-bielefeld.de), and 'Expiration' (31.08.2021). Below this is the 'Mass Spec Submission Form' and 'NMR Submission Form' tabs. The 'Sample Data' section includes 'Sample Name' (MP-01), 'Purification' (- none - / mixture / raw), and 'melts/boils' (don't know). There are checkboxes for 'Sample is dissolved' and 'Empirical Formula is unknown!'. The 'Empirical Formula' field contains 'C8H10N4O2' and is highlighted with a green box labeled 'a)'. To the right, 'Monoisotopic Mass (LUPAC)' and 'Monoisotopic Mass (lowest)' are both listed as 194,0804. The 'Toxicity (4)' section has checkboxes for 'carcinogenic', 'explosive', 'irritant', 'toxic', and 'untested' (checked). The 'Properties (4)' section has checkboxes for 'oxidation sensitive', 'small amount (<100µg)', 'strong odor', and 'temperature sensitiv'. The 'Sample is soluble in: (4)' section has checkboxes for 'Dichloromethane (Cl2)', 'Ethanol (EtOH)', 'Ethylacetate (EtOAc)', 'Isopropanol (IsoPrOp)', 'Methanol (MeOH)', and 'Tetrahydrofuran (THF)'. The 'Analytical Requirements' section has a 'Request for: (2)' list with 'Absence of Compound', 'Accurate Mass (±5 ppm)' (checked), 'Fragmentation Pattern', and 'MW (±0.2 Da)'. This list is highlighted with an orange box labeled 'b)'. To the right, 'Analogous samples have successfully been measured with: (2)' includes 'CI (Chemical Ionization)', 'EI (Electron Ionization)', 'ESI (Electrospray Ionization)' (checked), and 'MALDI (Matr. Ass. Laser Des. Ion.)'. The 'Chromatography necessary: (1)' section has checkboxes for 'GC-MS' (checked), 'HPLC-MS', and 'TLC Extraction'. Below these are three informational boxes: 'Request Info - Please consider!', 'Method Info - Please consider!' (with text: 'ESI: suited for polar, high boiling compounds, even salts'), and 'Separation Info - Please consider!'. A 'Comment/Info for the Operator' field contains 'Pharmakologisch aktiv!'. At the bottom, there is a 'Small Reaction Scheme Preview' showing a chemical structure, a 'Create MS Sample Submission Form' button, and a warning: 'Do not resize in drawing software before insert (check quality with preview!)'. There are also buttons for 'Insert Reaction Scheme from Clipboard' and 'Clear'.

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Mass Spec Submission Form NMR Submission Form

Sample Data

Sample Name: MP-01 Purification: -none - / mixture / raw melts/boils: don't know

Sample is dissolved

Empirical Formula is unknown! a)

Approx. Molecular Weight: 335.2 Da

Toxicity (4): carcinogenic, explosive, harmful, irritant, toxic, untested

Properties (4): oxidation sensitive, small amount (<100g), stable, strong odor, temperature sensitive

Sample is soluble in: (4): Dichloromethane (Cl), Ethanol (EtOH), Ethylacetate (EtOAc), Isopropanol (IsoProp), Methanol (MeOH), Tetrahydrofuran (THF)

Analytical Requirements

Request for: (2): Absence of Compound, Accurate Mass (+/- 5 ppm) b), Fragmentation Pattern, MW (+/- 0.2 Da)

Analogue samples have successfully been measured with: (2): CI (Chemical Ionization), EI (Electron Ionization), ESI (Electrospray Ionization), MALDI (Matr. Ass. Laser Des. Ion.)

Chromatography necessary: (1): GC-MS, HPLC-MS, TLC Extraction

Request Info - Please consider! Method Info - Please consider! Separation Info - Please consider!

Routine: EST: suited for polar, high boiling compounds, even salts

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Solvens	ESI	MALDI	EI
Acetonitril	●	●	●
Wasser (<80%!)	●	●	●
Methanol	●	●	●
Isopropanol	●	●	●
Tetrahydrofuran	●	●	●
Aceton	●	●	●
Chloroform	●	●	●
Dichlormethan	●	●	●
Toluol	●	●	●
Cyclohexan	●	●	●
DMSO	●	●	●
DMF	●	●	●
Pyridin	●	●	●

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Tutorial – MS request generation with SpecQuest

Table of Contents

Table of Contents	9
Tutorial – MS requests using SpecQuest	10
General remarks	10
1. Standard measurement (MW +/- 0.2 Da) – anticipated compound (e.g. product of synthesis)	11
2. Standard measurement (MW +/- 0.2 Da) – unknown compound.....	12
3. Accurate mass determination (MW +/- 5 ppm) – anticipated compound 13	
4. Accurate Mass determination (MW +/- 5 ppm) – unknown compound..	14
Table 1: Overview of compatible solvents for MS analytics	15

Tutorial – MS requests using SpecQuest

General remarks:

- I.) SpecQuest-Software can be downloaded from the [homepage of the MS department](#).
- II.) Requests for new SpecQuest Accounts or for renewal of existing ones can be made by email to Dr. Jens Sproß (name, email address, workgroup, phone number lab, duration of employment).
- III.) Please supply your sample as a solid. Sample solutions are only accepted in exceptional cases after consultation with the staff members of the MS department.
- IV.) Samples must not contain salts and buffers. These hinder the analysis (signal suppression) and can render the analysis impossible. Additionally, there is the possibility that the mass spectrometer will become unavailable due to contamination. Even residual amounts of hardly volatile solvents (DMF, DMSO) or other substances (e.g. TFA in case of ESI-measurements) can hamper the analysis.
- V.) In case of very sensitive samples, please consult the staff members of the MS department for an appointment to deliver and analyze your samples. If dry solvents are recommended for the analysis of your sample, please supply a sufficient amount (~5 mL/sample).
- VI.) Consult [table 1, page 15](#) for further advice regarding MS-compatible solvents. **Samples dissolved in DMSO, DMF or Pyridine, will not be processed!**
- VII.) Samples are submitted in the chemical cabinet on F02 next to room 217. Please put your MS-request on the corresponding shelf.
- VIII.) Return of your sample and of the analysis results is also performed via the chemical cabinet on F02 next to room 217. You are required to take your sample containers with you after the analysis is performed.
- IX.) MS results are supplied as print-out. Additionally, spectra of most measurements are available on the spectra-server of the faculty for chemistry.
- X.) In case of questions to MS-experiments and special measurements (e.g. proteins), please consult the staff members of the MS department.
- XI.) For citation of MS results, please consider the information in this [document](#). In case of questions, please consult Dr. Jens Sproß.

Contact persons:

Dr. Jens Sproß – all measurements, protein analytics, data evaluation, citation of MS methods and MS results in publications.

Sandra Heitkamp – measurements using the Esquire3000.

Ute Koppenbrink - Messungen am Esquire3000 and Synapt.

1. Standard measurement (MW +/- 0.2 Da) – anticipated compound (e.g. product of synthesis)

Creation of a MS request for a product of synthesis or compound with known sum formula.

- a) Switching between MS request and NMR request.
- b) Select your work group and name. **Please notice:** the date of expiry of your SpecQuest account is given in the area right to the selection field. Please ask for a renewal of your account in due time (send email to Dr. Jens Sproß).
- c) Fill in your sample name. Do not use a too long sample name (max. 10 characters). **Please do not use special characters and dots, as this will lead to problems during (automated) data acquisition.**
- d) Give further information regarding your sample (e.g. sample clean-up).
- e) Specify the sum formula of your compound (neutral compound. Please do not add charge carriers!). In case of ionic compounds (e.g. complexes), please indicate only the sum formula of the relevant ion, e.g. the cation.
- f) Indicate the solvent, which is able to dissolve your compound. **Please notice:** A table of the usual solvents that are compatible with the respective MS method can be found on [page 15](#).
- g) Select the type of MS-experiment. The standard measurement is MW +/- 0.2 Da. Please select “Fragmentation pattern” or “absence of compound” only when necessary. Consult the staff members of the MS department in such cases prior to your sample submission. For MS requests for “Accurate Mass”, please refer to [section 3, page 13](#).
- h) Indicate an ionization method, which was employed successfully in previous measurements of similar compounds.
- i) Please use the “comment for operator” for additional hints/informations (e.g. impurities, further relevant analytes in difficult to separate mixtures, counterion, ...).
- j) Insert a structure of your compound this can be achieved by copy/paste from common molecule editors. Please check the print-out for sufficient quality and readability. **Please notice:** The structure of a compound provides additional information and helps us to conduct the analysis.
- k) Print your MS request, remove the name tag for your sample (bottom right) and attach it to the sample container.

SpecQuest

Exit Server Info

Group OC1 Phone 1234
Name Person, Muster e-Mail muster.person@uni-bielefeld.de
Expiration 31.08.2021

Mass Spec Submission Form NMR Submission Form

Sample Data

Sample Name MP-01 Purification - none - / mixture / raw melts/boils don't know

Sample is dissolved
 Empirical Formula is unknown!

Empirical Formula C8H10N4O2 Monoisotopic Mass (IUPAC) 194,0804 Monoisotopic Mass (lowest) 194,0804

Toxicity (4) Properties (4) Sample is soluble in: (4)

carcinogenic oxidation sensitive
 explosive small amount (<100µg)
 irritant strong odor
 toxic temperature sensitive
 untested

Analytical Requirements

Request for: (2) Analogous samples have successfully been measured with: (2) Chromatography necessary: (1)

Absence of Compound CI (Chemical Ionization)
 Accurate Mass (±5 ppm) EI (Electron Ionization)
 Fragmentation Pattern ESI (Electrospray Ionization)
 MW (±0.2 Da) MALDI (Matr. Ass. Laser Des. Ion.)

Request Info - Please consider! Method Info - Please consider! Separation Info - Please consider!

Provide previously meas. spectrum! Sample must be chem. pure!

ESI: suited for polar, high boiling compounds, even salts

Comment/Info for the Operator
pharmacologically activ!

Small Reaction Scheme Preview

Create MS Sample Submission Form

Do not resize in drawing software before insert!
(check quality with preview!)
Insert Reaction Scheme from Clipboard
Clear

2. Standard measurement (MW +/- 0.2 Da) – unknown compound

Creation of a MS request for a product of synthesis or compound with unknown sum formula.

The creation of a MS request for the analysis of a compound with unknown sum formula is very similar to section 1. Only the differences are discussed here. Please proceed accordingly to [section 1, page 11](#) for the remaining part.

- Select "unkown compound".
- Indicate the approximate mass of your compound as accurate as possible (e.g. 400-500 Da). Please avoid to large mass ranges possible (e.g. 100-1000 Da).
- For further notes/informations (e.g. substance class) use this field.
- Otherwise, proceed in the same way as described in [section 1, page 11](#). For the determination of accurate masses of unknown compounds, please refer to [section 4, page 14](#).

SpecQuest

Exit Server Info

Group: OC1 Phone: 1234
Name: Person, Muster e-Mail: muster.person@uni-bielefeld.de
Expiration: 31.08.2021

Mass Spec Submission Form NMR Submission Form

Sample Data

Sample Name: MP-01 Purification: - none - / mixture / raw melts/boils: don't know

Sample is dissolved

Empirical Formula is unknown! **a)**
Approx. Molecular Weight: 300-350

Toxicity (4):
 carcinogenic
 explosive
 harmful
 irritant
 toxic
 untested

Properties (4):
 oxidation sensitive
 small amount (<100r
 stable
 strong odor
 temperature sensitiv

Sample is soluble in: (4)
 Dichloromethane (Cl
 Ethanol (EtOH)
 Ethylacetate (EtOAc
 Isopropanol (IsoProp
 Methanol (MeOH)
 Tetrahydrofuran (TH

Analytical Requirements

Request for: (2) **b)**
 Absence of Compound
 Accurate Mass (±5 ppm)
 Fragmentation Pattern
 MW (±0.2 Da)

Analogous samples have successfully been measured with: (2)
 CI (Chemical Ionization)
 EI (Electron Ionization)
 ESI (Electrospray Ionization)
 MALDI (Matr. Ass. Laser Des. Ion.)

Chromatography necessary (1)
 GC-MS
 HPLC-MS
 TLC Extraction

Request Info - Please consider! Method Info - Please consider! Separation Info - Please consider!

Routine: ESI: suited for polar, high boiling compounds, even salts

Comment/Info for the Operator **c)**
NMR gives hints to approx. 15-20 C-atoms. Contains probably O and N.

Small Reaction Scheme Preview

Create MS Sample Submission Form

Do not resize in drawing software before insert (check quality with preview!)
Insert Reaction Scheme from Clipboard
Clear

3. Accurate mass determination (MW +/- 5 ppm) – anticipated compound

Creation of a MS request for a product of synthesis or compound with known sum formula.

Important: An exact mass can only be determined after a standard measurement has been carried out (see [section 1, page 11](#)).

If you know in advance that you need an accurate mass determination, please *submit two MS requests (standard measurement **and** exact mass)* for your substance. If the standard measurement is successful, we carry out the determination of the accurate mass. *This saves you the hassle of bringing measurement requests and samples to F02 twice.* The MS department saves solvents and time for sample preparation, and there is also a lower risk potential, since the sample preparation does not have to be carried out twice!

If your substance has already been analyzed before, please enclose the complete measurement results with the application for accurate mass determination.

Please only use paper clips to combine several documents. Do not staple measurement requests and/or measurement results together (see note on the printed MS requests).

a) Specify the sum formula of your compound (neutral compound. Please do not add charge carriers!). In case of ionic compounds (e.g. complexes), please indicate only the sum formula of the relevant ion, e.g. the cation.

b) Select “Accurate Mass, +/- 5 ppm” as MS experiment.

c) Otherwise, proceed in the same way as described in [section 1, page 11](#).

Please notice: only one name tag has to be fixed to the sample container.

SpecQuest

Exit Server Info

Group: OC1 Phone: 1234
Name: Person, Muster e-Mail: muster.person@uni-bielefeld.de
Expiration: 31.08.2021

Mass Spec Submission Form NMR Submission Form

Sample Data

Sample Name: MP-01 Purification: - none - / mixture / raw melts/boils: don't know

Sample is dissolved

Empirical Formula is unknown! a)
Empirical Formula: C8H10N4O2
Monoisotopic Mass (IUPAC): 194,0804
Monoisotopic Mass (lowest): 194,0804

Toxicity (4):
 carcinogenic
 explosive
 irritant
 toxic
 untested

Properties (4):
 oxidation sensitive
 small amount (<100r
 stable
 strong odor
 temperature sensitive

Sample is soluble in: (4)
 Dichloromethane (Cl)
 Ethanol (EtOH)
 Ethylacetate (EtOAc)
 Isopropanol (IsoProp)
 Methanol (MeOH)
 Tetrahydrofuran (TH)

Analytical Requirements

Request for: (2) b)
 Absence of Compound
 Accurate Mass (+/- 5 ppm)
 Fragmentation Pattern
 MW (±0.2 Da)

Analogous samples have successfully been measured with: (2)
 CI (Chemical Ionization)
 EI (Electron Ionization)
 ESI (Electrospray Ionization)
 MALDI (Matr. Ass. Laser Des. Ion.)

Chromatography necessary: (1)
 GC-MS
 HPLC-MS
 TLC Extraction

Request Info - Please consider!
Provide previously meas. spectrum! Sample must be chem. pure!

Method Info - Please consider!
ESI: suited for polar, high boiling compounds, even salts

Separation Info - Please consider!

Comment/Info for the Operator
pharmacologically active!

Small Reaction Scheme Preview
Create MS Sample Submission Form

Do not resize in drawing software before insert (check quality with preview)!

Insert Reaction Scheme from Clipboard
Clear

4. Accurate Mass determination (MW +/- 5 ppm) – unknown compound

Creation of a MS request for a product of synthesis or compound with unknown sum formula.

Important: An exact mass can only be determined after a standard measurement has been carried out (see [section 2, page 12](#)).

Please enclose the complete measurement results to the application for accurate mass determination. The accurate mass cannot be determined without measurement results from the standard measurement!

Please only use paper clips to combine several documents. Do not staple measurement requests and/or measurement results together (see note on the printed MS requests).

- Indicate the mass of your compound. **Important:** mass of the neutral compound, *not as adduct with an ion, no mass range!*
- Select “Accurate Mass, +/- 5 ppm” as MS experiment.
- Please specify uncommon elements in your compound. For the generation of the analysis result only the elements C, H, N, and O are considered, if not stated otherwise. Depending on the adduct-ion-species, Na and Cl are additionally taken into account.
- Otherwise, proceed in the same way as described in [section 1, page 11](#).

The screenshot shows the SpecQuest software interface. At the top, there are fields for Group (OC1), Name (Person, Muster), and Phone (1234). Below this is the 'Mass Spec Submission Form' section. The 'Sample Data' section includes 'Sample Name' (MP-01), 'Purification' (- none - / mixture / raw), and 'melts/boils' (don't know). A green box highlights the 'Empirical Formula is unknown!' checkbox and the 'Approx. Molecular Weight' field (335.2). The 'Analytical Requirements' section has a 'Request for:' section with 'Accurate Mass (+/- 5 ppm)' selected, highlighted by an orange box. The 'Comment/Info for the Operator' field at the bottom contains the text 'NMR gives hints to approx. 15-20 C-atoms. Contains probably O and N.', highlighted by a red box. At the bottom of the interface, there is a 'Small Reaction Scheme Preview' showing a chemical structure and a button to 'Insert Reaction Scheme from Clipboard'.

You will receive a table of potential sum formulas. We will be happy to help you filter this list for the most likely sum formula. Please contact Dr. Jens Sproß about this.

Table 1: Overview of compatible solvents for MS analytics

Solvent	ESI	MALDI	EI
Acetonitrile	●	●	●
Water (<80%!)	●	●	●
Methanol	●	●	●
Isopropanol	●	●	●
Tetrahydrofuran	●	●	●
Acetone	●	●	●
Chloroform	●	●	●
Dichlormethane	●	●	●
Toluene	●	●	●
Cyclohexane	●	●	●
DMSO	●	●	●
DMF	●	●	●
Pyridine	●	●	●

Key:


- - compatible
- - partially possible, e.g. as part of a solvent mixture in ESI and MALDI
- - not compatible to this method

General note: Samples dissolved in DMSO, DMF and pyridine will not be processed! Please remove these solvents before submitting your samples.

Note ESI: In case of ESI measurements, adding a small amount of formic acid, NaClO₄, LiCl, or AgBF₄, respectively, the analyte can be ionized better. This will lead to the formation of H⁺-, Na⁺-, Li⁺- or Ag⁺-adduct-ions, respectively. The addition is noted on the print-outs of the spectra. Rule of thumb: addition of formic acid supports the ionization of nitrogen-containing analytes, addition of NaClO₄ supports the ionization of oxygen-containing analytes and addition of AgBF₄ supports the ionization of aromatic and polyunsaturated hydrocarbons.

IMPORTANT: *These substances are only added by the employees of the mass spectrometry department as part of the sample preparation!*

Note MALDI: for MALDI measurements, different matrices are used depending on the analyte: CHCA (α -cyano-4-hydroxycinnamic acid) for small molecules and peptides, DHB (2,5-dihydroxybenzoic acid) for small molecules and proteins, and DCTB (*trans*-2-[3-(4-*tert*-butylphenyl)-2-methyl-2-propenylidene]-malonodinitril) for hydrocarbons. Rule of thumb: With CHCA and DHB, mainly H⁺ and Na⁺ adduct ions are detected, with DCTB radical cations!

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Mass Spec Submission Form | NMR Submission Form

Sample Data
Sample Name: MP-01 Purification: - none - / mixture / raw melts/boils: don't know

Sample is dissolved
 Empirical Formula is unknown! a)

Approx. Molecular Weight: 300-350 b)

Toxicity (4)
 carcinogenic
 explosive
 harmful
 irritant
 toxic
 untested

Properties (4)
 oxidation sensitive
 small amount (<100r
 stable
 strong odor
 temperature sensitiv

Sample is soluble in: (4)
 Dichloromethane (Cl
 Ethanol (EtOH)
 Ethylacetate (EtOAc)
 Isopropanol (IsoProp
 Methanol (MeOH)
 Tetrahydrofuran (TH

Analytical Requirements

Request for: (2)
 Absence of Compound
 Accurate Mass (± 5 ppm)
 Fragmentation Pattern
 MW (± 0.2 Da)

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Chromatography necessary (1)
 GC-MS
 HPLC-MS
 TLC Extraction

Request Info - Please consider
Routine

Method Info - Please consider
ESI: suited for polar, high boiling compounds, even salts

Separation Info - Please consider

Comment/Info for the Operator c)
NMR gives hints to approx. 15-20 C-atoms. Contains probably O and N.

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The screenshot shows the 'SpecQuest' software interface. At the top, there are fields for 'Group' (OC1), 'Name' (Person, Muster), 'Phone' (1234), 'e-Mail' (muster.person@uni-bielefeld.de), and 'Expiration' (31.08.2021). Below this is the 'Mass Spec Submission Form' section, which includes 'Sample Data' with 'Sample Name' (MP-01), 'Purification' (- none - / mixture / raw), and 'melts/boils' (don't know). The 'Empirical Formula' field is highlighted with a green box and contains 'CSH10N4O2', with a small 'a)' next to it. Below this are sections for 'Toxicity (4)', 'Properties (4)', and 'Sample is soluble in: (4)'. The 'Analytical Requirements' section is highlighted with an orange box and contains 'Request for: (2)' with 'Accurate Mass (±5 ppm)' selected, and 'b)' next to it. Other options include 'Absence of Compound', 'Fragmentation Pattern', and 'MW (±0.2 Da)'. The 'Method Info - Please consider!' section is also visible, with 'ESI: suited for polar, high boiling compounds, even salts' noted. At the bottom, there is a 'Small Reaction Scheme Preview' showing a chemical structure and a 'Create MS Sample Submission Form' button.

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Analytical Requirements

Request for: (2): Absence of Compound, Accurate Mass (+/- 5 ppm) b), Fragmentation Pattern, MW (+0.2 Da)

Analogous samples have successfully been measured with: (2): CI (Chemical Ionization), EI (Electron Ionization), ESI (Electrospray Ionization), MALDI (Matr. Ass. Laser Des. Ion.)

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Request Info - Please consider: Provide previously meas. spectrum! Sample must be chem. pure!

Method Info - Please consider: ESI: suited for polar, high boiling compounds, even salts

Separation Info - Please consider:

Comment/Info for the Operator c): NMR gives hints to approx. 15-20 C-atoms. Contains probably O and N.

Small Reaction Scheme Preview

Create MS Sample Submission Form

Do not resize in drawing software before insert (check quality with preview!)

Insert Reaction Scheme from Clipboard

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